



## Head Magnifier oem

[Read More](#)

**SKU:**

**Price:** 0.00 50

**Stock:** instock

**Categories:** [Devices](#)

### Product Description

2led headband illuminating magnifier with portable detachable led light source box and 11 kinds of magnifications



## Hygedent Dental Alginate Impression Material

[Read More](#)

**SKU:**

**Price:** 0.00 5

**Stock:** instock

**Categories:** [Consumables](#)

### Product Description

Dental Alginate Impression Material



## I-Gel Phosphoric Acid Etching Gel

[Read More](#)

**SKU:**

**Price:** 0.00 5

**Stock:** instock

**Categories:** [Consumables](#)

### Product Description

Phosphoric Acid Etching Gel



## IMD NITI FILES M-PRO

[Read More](#)

**SKU:**

**Price:** ₹.₹ 10

**Stock:** instock

**Categories:** [Consumables](#)

### Product Description

NITI FILES M-PRO SIZE : #18 #20 #25,25mm



## IMED Dental Cotton Roll

[Read More](#)

**SKU:**

**Price:** ₹.₹ 8

**Stock:** instock

**Categories:** [Disposables](#)

### Product Description

Cotton Roll SIZE 10×38MM Packing Quantity 2000PCS



## IMED Disposable Head Cap

[Read More](#)

**SKU:**

**Price:** ₹.₹ 3

**Stock:** instock

**Categories:** [Disposables](#)

### Product Description

Disposable Head Cap 100 PCS



## IMED Disposable Show Cover

[Read More](#)

**SKU:**

**Price:** ₺.₺ 3

**Stock:** instock

**Categories:** [Disposables](#)

### Product Description

Disposable Show Cover size 15×41 cm 100 PCS



## IMED Sterilization Rolls 75mm×200m

[Read More](#)

**SKU:**

**Price:** ₺.₺ 8

**Stock:** instock

**Categories:** [Disposables](#)

### Product Description

Sterilization Rolls 75mm×200m



## Implant Contra Angel 20 TO 1

[Read More](#)

**SKU:**

**Price:** ₺.₺ 120

**Stock:** instock

**Categories:** [Tools](#)

## Product Description

**Dental Implant Surgical handpiece 20:1 Reduction contra angle** Features: 1. 20:1 Reduction . For Ø2.35mm CA Bur 2. Driven by Implant Surgical Unit 3. Single Water Spray & Non-Optic 4. AI DENTAL stainless steel bearings 5.135°C Autoclavable 6. Warranty:one year for handle; 6 months for cartridge Performance specifications : connect: E-Type Implant Surgical motor Max. Rotation Speed (Motor) □40,000min-1 Max. Rotation Speed (Handpiece) □2,000 min-1 Max. Torque : 70 Ncm Gear Ratio □20:1 Reduction Bur Type□ Ø2.35mm CA Bur Optic □Non-Optic Coolant Water Type □Double Water Injection Packing List : Handpiece \*1 + Wrench \*1 + manual \*1 Packing Weight : 0.15kg Packing Size : 15\*7\*3 mm



## Interproximal reduction kit

IPR kit has all the specific tools for interproximal reduction.

[Read More](#)

**SKU:**

**Price:** ~~₹.₹200~~ Original price was: 200 ₹.₹.₹.₹ 150 Current price is: 150 ₹.₹.

**Stock:** instock

**Categories:** [Devices](#)

**Tags:** [inter proximal](#), [orthodontic](#), [reciprocation](#), [reduction](#)

## Product Description

This kit contains all necessary instruments used in interproximal reduction procedures.

## Contents And Features Of Interproximal Reduction Kit

IPR kit has all the specific tools for interproximal reduction. Interproximal reduction kit has a reciprocating handpiece which is used for mechanical seating of restorations during cementation. It also has strip stands with five holes for holding strips. Interproximal strip gauge that provides the right amount of stripping with labeled dimensions. And finally one chuck key. It has uniquely designed slow speed handpiece and other tools for IPR. All tools available in the IPR kit are autoclavable and easy to use.

## Uses And Indications of Interproximal reduction

Process of interproximal reduction is used for several procedures like creating space to improve tooth alignment. It helps to improve occlusion. This also helps to increase the stability of teeth for long term retaining of their position. IPR also helps to reduce black triangles between the teeth. It can be used to treat relapsed crowding. Basically interproximal reduction is used to treat malocclusions where extraction is not necessary. Black triangles which happen because of the absence of gingival papilla are also reduced by IPR. Stripping is advisable to realign anterior teeth without labial tipping. It is also used for cosmetic recontouring of the canines to mimic lateral incisors.

# Kidney Dish

[Read More](#)

**SKU:**

**Price:** ₹.₹ 8

**Stock:** instock

**Categories:** [Tools](#)



## Product Description

Kidney Dish



# Kit Probe And Tweezer

[Read More](#)

**SKU:**

**Price:** ₹.₹ 5

**Stock:** instock

**Categories:** [Tools](#)

## Product Description

Kit Probe And Tweezer