

Surgery Motor

[video width="852" height="480"

mp4="https://ecodentstores.com/wp-content/uploads/2024/09/Self-Pum ping-Cooling-Water-Dental-Implant-Brushless-Motor-Irrigation-Surgical.mp4"][/video]

Read More

SKU:

Price: [].[] 195
Stock: instock
Categories: Devices

Product Description

Dental Surgical and implant brushless motor system Brushless motor Torque: 17N.cm for 1:1 contra angle, and 50Ncm for 20:1 contra angel Rotation speed 0-50,000r/min Cooling water: 0-100ml/min Self-pumping water **Specification:**

Motor type Direct current brushless motor

Voltage DC36V Operation mode REV

Speed 2000-50,000 rpm

Torque 17N·cm

Sterilization 121°C20 minute, or 135°C12 minute

Motor size Φ23×L145mm

Noise ≤60dB

Control unit:

Model EM1-A Voltage 110-220V□

Motor size L225×W90×H100mm

Fuse protector T3.15AL250V

Water spray:

Cooling water supply 0-100ml/min

Maintenance: Use alcohol to wipe the motor and control panel everyday to keep it clean. It is absolutely forbidden to use oxidized ionized water (strongly oxidized water, peroxide water) or disinfectant for cleaning, soaking and wiping.