



## Air Propy Sander Gun

[Read More](#)

**SKU:**

**Price:** ~~30~~ Original price was: 30 ~~25~~ Current price is: 25

**Stock:** instock

**Categories:** [Devices](#)

### Product Description

Connection: 4 OR 2 hole Water Pressure: 0.05MPa - 0.5MPa Air Pressure: 0.2MPa-0.5MPa Max. Waterflow: 30ml/min Max. Powderflow: 2g-3g/min



## CURING LIGHT 5 SEC Cure

[Read More](#)

**SKU:**


**Price:** 25

**Stock:** instock

**Categories:** [Devices](#)

### Product Description

Medical Dental Curing Light 5 Seconds Cure



## Dental Acrylic Organizer

[Read More](#)

**SKU:**

**Price:** ~~30~~ Original price was: 30 ~~25~~ Current price is: 25

**Stock:** instock

**Categories:** [Devices](#)

## Product Description

Description 100% brand new and high quality High quality acrylic imported from UK. Transparent for placing resins, adhesives, syringes, applicators and other small tools. Net weight: 650g Packing List 1pc new Dental Acrylic Holder



## Dental led loupes

[Read More](#)

**SKU:**

**Price:** ~~£125~~ Original price was: 125 £. Current price is: 75 £.

**Stock:** instock

**Categories:** [Devices](#)

## Product Description

elegant dental loupes with led light . magnification: 2.5x distance :45 cm



## DTE O-Light Curing Light

[Read More](#)

**SKU:**

**Price:** ~~£90~~ Original price was: 90 £. Current price is: 85 £.

**Stock:** instock

**Categories:** [Devices](#)

## Product Description

The **DTE O-Light Curing Light** is a powerful tool used in dentistry for light-curing procedures. Here are the key details about this device:

- **Modes:**

- **High Mode:** Light intensity of **2300 mw/cm<sup>2</sup> - 2500 mw/cm<sup>2</sup>**. Suitable for orthodontics, porcelain veneers, and post-core restoration.
- **Standard Mode:** Light intensity of **1000 mw/cm<sup>2</sup> - 1200 mw/cm<sup>2</sup>**. Suitable for general curing.

- **Features:**

- **Fast Curing:** The DTE O-Light can cure within **1 second<sup>1</sup>**.
- **Maximum Power:** It solidifies **2 mm of composite per second**.
- **Ergonomic Design:** Lightweight and features a **360° head rotation**.
- **Wavelength:** Covers the range from **385 nm to 515 nm**, allowing it to polymerize all composites on the market.
- **Battery:** Rechargeable with a capacity of **2000 mAh** and overcharge protection.
- **Intensity:** Ranges from **1000 mw/cm<sup>2</sup> to 2500 mw/cm<sup>2</sup>**.
- **Transformer Voltage:** Compatible with **100V to 240V** input.
- **Metal Tip:** The LED is located at the end of the metal head, concentrating the light beam for effective curing.

- **Guarantee:** Comes with a **24-month warranty** (6 months for accessories).

The DTE O-Light is an essential tool for dental professionals, especially during procedures like bracket placement, porcelain veneers, and adhesives. Its rapid curing capability and precision make it a reliable choice for practitioners worldwide



## Endo Smart Endodontic Motor iq100

[video width="1280" height="720" mp4="https://ecodentstores.com/wp-content/uploads/2024/09/Dental-Wireless-Endo-Motor-Smart-with-LED-Lamp-16-1-Standard-Contra-Angle-Endodontic-Instrument-AliExpress-66-2.mp4"]

[Read More](#)

**SKU:**

**Price:** ₹.₹ 95

**Stock:** instock

**Categories:** [Devices](#)

### Product Description

#### 9 Different Programs:

- **Auto-reverse And Auto Forward Mode.**
- **The Slim Body and Extremely Small Head Allow Optimal Visibility During Procedures.**
- **Compatible with 16:1 dental handpiece**
- **Build-in Rechargeable Battery, over hours of continuous working**
- **Stable And Prolonged usage.**
- **High Brightness LED To Provide Better Illumination During RCT.**

#### TECHNICAL SPECIFICATIONS:

- Battery Voltage - DC 3.7V
- Battery Capacity - 950mAh
- Rated input Voltage-AC 110-220V 60/60HZ
- Preset Programs -9 Types
- Torque Range(N.cm)-0.3-3.0
- Speed Range(rpm)-140-550
- Working modes - 4 Kinds



## Gutta Cutter

[Read More](#)

**SKU:**

**Price:** ₹.₹ 20

**Stock:** instock

**Categories:** [Devices](#)

### Product Description

used for cutting the extra gutta percha from the root canal .

## Head Magnifier



[Read More](#)

**SKU:**

**Price:** ₹.₹ 25

**Stock:** instock

**Categories:** [Devices](#)

### Product Description

compact Head magnifier with lenses and led light

## Head Magnifier oem



[Read More](#)

**SKU:**

**Price:** ₹.₹ 50

**Stock:** instock

**Categories:** [Devices](#)

### Product Description

2led headband illuminating magnifier with portable detachable led light source box and 11 kinds of magnifications



## Interproximal reduction kit

IPR kit has all the specific tools for interproximal reduction.

[Read More](#)

**SKU:**

**Price:** ₹.₹-200 Original price was: 200 ₹.₹.₹.₹ 150 Current price is: 150 ₹.₹.

**Stock:** instock

**Categories:** [Devices](#)

**Tags:** [inter proximal](#), [orthodontic](#), [reciprocation](#), [reduction](#)

### Product Description

This kit contains all necessary instruments used in interproximal reduction procedures.

### Contents And Features Of Interproximal Reduction Kit

IPR kit has all the specific tools for interproximal reduction. Interproximal reduction kit has a reciprocating handpiece which is used for mechanical seating of restorations during cementation. It also has strip stands with five holes for holding strips. Interproximal strip gauge that provides the right amount of stripping with labeled dimensions. And finally one chuck key. It has uniquely designed slow speed handpiece and other tools for IPR. All tools available in the IPR kit are autoclavable and easy to use.

### Uses And Indications of Interproximal reduction

Process of interproximal reduction is used for several procedures like creating space to improve tooth alignment. It helps to improve occlusion. This also helps to increase the stability of teeth for long term retaining of their position. IPR also helps to reduce black triangles between the teeth. It can be used to treat relapsed crowding. Basically interproximal reduction is used to treat malocclusions where extraction is not necessary. Black triangles which happen because of the absence of gingival papilla are also reduced by IPR. Stripping is advisable to realign anterior teeth without labial tipping. It is also used for cosmetic recontouring of the canines to mimic lateral incisors.



## Marathon Micromotor

[Read More](#)

**SKU:**

**Price:** ₹.₹ 120

**Stock:** instock

**Categories:** [Devices](#)

### Product Description

Marathon Micromotor



# Micropex Apex Locator

[Read More](#)

**SKU:**

**Price:** ₱.₱ 50

**Stock:** instock

**Categories:** [Devices](#)

## Product Description

Measuring Wire File Clip Lip Hook Touch Probe Tester Battery : 3pcs of AAA Battery Power : 4.5V 0.1A